

WARNINGS

Before performing one or more treatments with the CTU Mega 20 Diamagnetic Pump device, you must submit the patient to a medical screening in order to detect:

- Unsuitable physiological states
- Presence of ferromagnetic material within the areas of the body to be treated

CAUTION: Ensure that both the operator and the patient do not wear any metallic object.

The application is absolutely prohibited on patients/subjects with:

- Serious Heart and kidney disease
- Fertile epiphysis in the treatment area
- Pregnancy (fetus in the treatment area)
- Tumors in the treatment area and malignancies

For other contraindications, precautions and instructions for use of the device read the supplied User Manual.